

LAST WEEK IN THE POTIONS LABORATORY

Separated mextures based on deffering physical and chemical properties

Used Excel to calculate average, standard deveation, and weight percents

Q≒test in action!

BUCHNER FİLTERİNG

EXPERÎMENT 5
QUALÎTATÎVE CHEMÎCAL ANALYSÎS

PREVYOUS YEARS

Gorgon's blood Iặquid gold Vitreous humor of a blind mule



Ac†ds Bases Inorgan†c salt

"QUANT" VS "QUAL"

QUANTÉTATÉVE – HOW MUCH ÉS THERE?

QUALÉTATÉVE – WHAT ÉS ÉT?

You will identify the chemical identities of 5 unknown solutions based on how they react (or don't react!)

with one another.

TWO-PART LAB

- PART 1: Mix ten known solutions and record the results of the deactions
- PART 2: Mix your five unknowns and compare the results with what you saw in Part One.

BE EXACT!	
THE MORE ACCURATELY YOU RECORD YOUR OBSERVATIONS, THE EASTER IT WILL BE TO IDENTIFY YOUR UNKNOWNS.	
H . T-	
The Known Solutions	
ACTOS: HCl H ₂ SO ₄ HNO ₃	
BASES: NaOH Na ₂ S	
SALTS: Ba(NO_3) ₂ Ag NO_3 K ₂ CrO ₄	
$Fe(NO_3)_3 Ni(NO_3)_2$	
Ave coverage and O 10 on O 20 A	
ALL SOLUTIONS ARE Q.1Q OR Q.2Q M.	
What are we looking for?	

*Explosions	
*Supernatural creatures	
*Rashes, mutations,	
TRANSFORMATŤONS	

What are we'llooking for?

PRECÉPÉTATES, MOSTLY.
[SEE THE SOLUBÉLÉTY TABLE ÉN APPENDÉX 2.]

Don't expect to see any actombase action.

RŤDDLE ME THŤS

What do you get when you cross hydrochlorec aced with selver netrate?

Balanced Chemical Equation

$$HCl + AgNO_3 \longrightarrow AgCl + HNO_3$$

ADD THE PHYSICAL STATES OF EACH COMPOUND

$$HCl_{(aq)} + AgNO_{3^{(aq)}} \longrightarrow AgCl_{(s)} + HNO_{3^{(aq)}}$$

This is sometimes called a molecular equation.

LET'S GET REAL

$$HCl_{(aq)} + AgNO_{3^{(aq)}} \longrightarrow AgCl_{(s)} + HNO_{3^{(aq)}}$$

Total Ionec Equateon – wrete aqueous compounds as endevedual eons:

LOTS OF SPECTATOR FONS

Time to Clean House

Cross out spectator fons

$$H^{+}_{(aq)} + Cl^{-}_{(aq)} + Ag^{+}_{(aq)} + NQ_{3}^{-}_{(aq)} \longrightarrow AgCl_{(s)} + H^{+}_{(aq)} + NQ_{3}^{-}_{(aq)}$$

NET IONEC EQUATEON

$$Ag^{+}_{(aq)} + Cl^{-}_{(aq)} \longrightarrow AgCl_{(s)}$$

THE NET IONEC EQUATEON

$$NaCl_{(aq)} + AgNO_{3^{(aq)}} \longrightarrow NaNO_{3^{(aq)}} + AgCl_{(s)}$$

$$Ba(Cl)_{2} \stackrel{\text{(aq)}}{=} + 2AgCH_{3}COO_{\text{(aq)}} \longrightarrow Ba(CH_{3}COO)_{2} \stackrel{\text{(aq)}}{=} + 2AgCl_{\text{(s)}}$$

$$NH_4Cl_{(aq)} + AgClO_{3^{(aq)}} \longrightarrow NH_4ClO_{3^{(aq)}} + AgCl_{(s)}$$

All of these reactions have the same net ionic equation:

$$Ag^{+}_{(aq)} + Cl^{-}_{(aq)} \longrightarrow AgCl_{(s)}$$

STMPLE TS GOOD

- The net ionic fourtion describes the chemical reaction that occurs, and does not include any ions that do not take part in the reaction, even though those ions are present in solution.
- How do we know which fons will react and which ones won't?

SOME QUÝCK SOLUBĚLĚTY RULES

* All compounds containing alkali metals and ammonium ion are soluble.

 $Li^{\scriptscriptstyle +} \quad Na^{\scriptscriptstyle +} \quad K^{\scriptscriptstyle +} \quad Rb^{\scriptscriptstyle +} \quad Cs^{\scriptscriptstyle +} \quad NH_4^{\; \scriptscriptstyle +}$

* ALL COMPOUNDS CONTAÎNÎNG NÎTRATE, CHLORATE, DERCHLORATE, AND ACETATE ADE SOLURIE

 $NO_3^ \overline{ClO_3^-}$ $ClO_4^ \overline{CH_3COO^-}$ or $\overline{C_3H_3O_2^-}$

SOME QUÍCK INSOLUBILITY RULES

- $^{\circ}$ All compounds containing PO₄ 3 , CO $_3^2$ -, and SO $_3^2$ ions are insoluble except those that also contain alkali metals or NH₄ $^+$.
- * All hydroxides and sulfides are insoluble except those that also contain alkali metals, or NH₄⁺. Hydroxides end with –OH, and sulfides end with –S. (Exception: some group II hydroxides are slightly soluble.)
- When in doubt, $Ag^+ Pb^{2+}$ and Hg compounds tend to be insoluble.

IN THE POTIONS LABORATORY

- * Create an array of reactions in the microwell plate similar to the one in the lab manual.
- * Use only 2 = 5 drops of each reactant.
- * DO NOT touch the type of the dropper bottles to the solutions in the Microwell plate or you will die a most painful death.

Vile, hideous fluids! CHROMŤUM (VI) ŤN PARTŤCULAR ŤS NASTY! EMPTY YOUR USED MECROWELL PLATES THE DESCUSTING PLASTIC TRAY IN THE HOOD. PINSE THE PLATES WITH ORDINARY TAP WATER, DO A FÎNAL PÎNSE WÎTH DEÎON WATER, AND STACK THE PLATES IN THE HOOD. LAB REPORT MOLECULAR EQUATIONS FOR 11 PRECEPETATEON REACTEONS. NET IONEC FOUATEONS FOR 11 PRECEPETATEON REACTEONS. 11 * 11 is 22 fountions altogether. Post Lab Preview CaBr₂ (aq) **and** NaBr (aq) What will precipitate with Ca^{2+} BUT NOT WITH Na+? $Fe(NO_3)_3 \ {}_{\text{(aq)}} \ \textbf{AND} \ Mg(ClO_3)_2 \ {}_{\text{(aq)}}$ What will precipitate with Fe^{3+} but not with Mg^{2+} ?

NEXT WEEK

EXPERÎMENT 4: ACÎD BASE TÎTRATÎON

- MORE TOVOLVED THAN THIS LAB.
- Requires more preparation than our previous lars.
- START PREPARING FARIY

Quiz Time

- * You will need a calculator every week (except next week).
- * Make sure you know your section number and your TA's name.

6