


## Experiment 4 Acid-Base Titration

CH 204 Fall 2006  
Dr. Brian Anderson



---

---

---

---

---

---

---

---



### What We Lernd in Skool Last Week

Molecular Equations

Simple Solubility Rules

Spectator Ions and Net Ionic Equations

Microscale Techniques

---

---

---

---

---

---


---

---



### This Week: Acid-Base Titrations

What exactly *is* a titration, anyway?



At the equivalence point  
Moles  $\text{H}^+$  = Moles  $\text{OH}^-$

---

---

---

---

---

---

---

---

## Titration Setup

Burette containing **NaOH**.

Read all volumes to **0.01 mL!**

Erlenmeyer flask containing **acid sample, water, and two drops of phenolphthalein**.

---

---

---

---

---

---

---

---

## Phenolphthalein

A weak acid.

Colorless below pH 8, pink above pH 8.

Your acid solution will go from colorless to faint pink at the endpoint. If it turns bright pink, you have gone too far.

<http://www.chemistry.wustl.edu/~courses/genchem/Labs/AcidBase/phph.htm>

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---



## Review of pH

1

Neutral                    pH = 7  
 Acidic                    pH < 7  
 Basic (alkaline)        pH > 7

pH is like a Richter scale for acid concentrations, but the higher the pH, the lower the [H<sup>+</sup>].

---

---

---

---

---

---

---

---



## pH = log<sub>10</sub>[H<sup>+</sup>]

2

[H <sup>+</sup> ]	log[H <sup>+</sup> ]	pH
10 <sup>-1</sup>	-1	1
10 <sup>-2</sup>	-2	2
10 <sup>-3</sup>	-3	3
10 <sup>-4</sup>	-4	4
1.78 × 10 <sup>-5</sup>	-4.75	4.75
4.68 × 10 <sup>-7</sup>	-6.33	6.33
8.13 × 10 <sup>-12</sup>	-11.90	11.90

---

---

---

---

---

---

---

---



## Calculating pH

3

**Given the [H<sup>+</sup>], pH = -log[H<sup>+</sup>]**

What is the pH of a 0.025 M HCl solution?

$$\text{pH} = -\log[.025] = -(-1.6) = 1.6$$

**Given the pH, [H<sup>+</sup>] = 10<sup>-pH</sup>**

What is the [H<sup>+</sup>] concentration in human arterial blood, which has a pH of 7.40?

$$[\text{H}^+] = 10^{-\text{pH}} = 10^{-7.40} = 3.98 \times 10^{-8} \text{ M}$$

---

---

---

---

---

---

---

---



## Review

4

If you are given the  $[H^+]$

$$pH = -\log[H^+]$$

If you are given the pH

$$[H^+] = 10^{-pH}$$

---

---

---

---

---

---


---

---



## Today: Titration Marathon!

Determine the concentration  
of an unknown acidic solution:



$$\text{Moles } H^+ = \text{Moles } OH^-$$

$$M_{H^+} \times V_{H^+} = M_{OH^-} \times V_{OH^-}$$

$$M_{acid} \times V_{acid} = M_{base} \times V_{base}$$

---

---

---

---

---

---

---

---



## Experiment 4 Overview

### PART 1: STANDARDIZATION OF NaOH

Mix up 1 liter of NaOH solution.

Weigh out 2 grams of KHP powder,  
dissolve in 75 ml water, **ADD**  
**PHENOLPHTHALEIN**, and titrate (3×).

Calculate concentration of NaOH using

$$\text{Moles of Acid} = \text{Moles of Base}$$

---

---

---

---

---

---

---

---



## Moles solid = Moles aqueous

$$\frac{\text{Mass of KHP}}{\text{MW of KHP}} = M_{\text{NaOH}} \times V_{\text{NaOH}}$$

$$\frac{\text{Mass of KHP}}{\text{MW of KHP} \times V_{\text{NaOH}}} = M_{\text{NaOH}} \text{ (0.xxxx M)}$$

---

---

---

---

---

---



---

---



## KHPthalate

Molecular weight is NOT K + H + P = 71.




---

---

---

---

---


---

---

---



## KHPthalate vs Cthulhu



**KHPHTHALATE**



**CTHULHU**

---

---

---

---

---

---

---

---



## Part Two: A Return to the Potions Lab

Identify your unknown acid sample using the qualitative reactions from last week.



---

---

---

---

---

---

---

---



## Part 3: Titrate Your Unknown

5 ml unknown acid, 75 ml water, and 2 drops of phenolphthalein in a 250 ml flask.  
Titrate using NaOH (3×)

In an ideal world, you will get the exact same  $V_{\text{NaOH}}$  all three times.

Calculate the molarity of your acid.

---

---

---

---

---

---

---

---



## Moles<sub>H<sup>+</sup></sub> = Moles<sub>OH<sup>-</sup></sub>

For HCl and HNO<sub>3</sub>,

$$M_{\text{acid}} \times V_{\text{acid}} = M_{\text{base}} \times V_{\text{base}}$$

For H<sub>2</sub>SO<sub>4</sub>

$$2 \times M_{\text{acid}} \times V_{\text{acid}} = M_{\text{base}} \times V_{\text{base}}$$

$$V_{\text{acid}} = 5.00 \text{ ml}$$

---

---

---

---

---

---

---

---



## General Form for Acid-Base Titrations

Moles H<sup>+</sup> = Moles OH<sup>-</sup>

$$M_{\text{acid}} \times V_{\text{acid}} \times \# \text{ of H}^+ = M_{\text{base}} \times V_{\text{base}} \times \# \text{ of OH}^-$$

---

---

---

---

---

---

---

---



## Part 4: Citric Acid in Juice

Orange, Grapefruit, or Pineapple  
15 ml juice, 60 ml water, and  
**2 drops of phenolphthalein.**

Titrate just once. Endpoint is hard to see  
in orange juice.

---

---

---

---

---

---

---

---



## A word about citric acid

That word is triprotic!

**1** Mole of citric acid = **3** moles of H<sup>+</sup>

So the number of moles of H<sup>+</sup> is **3 times**  
the number of moles of citric acid:

$$\underline{3} \times M_{\text{Citric acid}} \times V_{\text{Citric acid}} = M_{\text{base}} \times V_{\text{base}}$$

---

---

---

---

---

---

---

---



**All your base are belong to us**

Leftover NaOH goes into the waste container in the hood.

Keep your unknown acid for now.

**DO YOUR CALCULATIONS BEFORE YOU DUMP YOUR LEFTOVER BASE!!**

If you have time, fill in all the data tables before you leave the lab.

---

---

---

---

---

---

---

---