

Name Mark PriorDate 09/14/05TA Carlos ZambranoSection # 52420**Title: "Are the Densities of Coke and Diet Coke Different?"**

Objective: To measure and compare densities of Coke and Diet Coke. To evaluate precision of the volume measurements performed with a 5 ml volumetric pipette, a 10 ml graduated cylinder, and a 50 ml burette. To verify that density is an intensive property.

Sample: Diet Coke

Sample temperature: 22.5 °C

Part 1. "Are the Densities of Coke and Diet Coke Different?"**My data**

| | Volumetric Pipette | | Graduated Cylinder | | Burette | |
|------------------------------------|--------------------|---------|--------------------|---------|---------|---------|
| | Trial 1 | Trial 2 | Trial 1 | Trial 2 | Trial 1 | Trial 2 |
| Mass of Empty Beaker (g) | 28.4576 | 36.1582 | 32.0014 | 37.2285 | 46.2828 | 30.6521 |
| Mass of Beaker + Liquid Sample (g) | 33.4530 | 41.1399 | 36.6603 | 42.0284 | 51.1820 | 35.5737 |
| Mass of Liquid Sample (g) | 4.9954 | 4.9817 | 4.6589 | 4.7999 | 4.8992 | 4.9216 |
| Volume Delivered (ml) | 5 | 5 | 5 | 5 | 5 | 5 |
| Calculated Density (g/ml) | 0.9991 | 0.9963 | 0.9318 | 0.9600 | 0.9798 | 0.9843 |

(NOTE: Numbers listed in the Table above should be used as an example ONLY. They do not represent a real set of data.)

Sample Calculation for Trial 1:

Mass of liquid sample = 33.4530 g - 28.4576 g = 4.9954 g

$$\text{Density} = \frac{m}{V} = \frac{4.9954 \text{ g}}{5 \text{ ml}} = 0.9991 \text{ g/ml}$$

Name Mark PriorDate 09/14/05TA Carlos ZambranoSection # 52420**Results from the combined class data**

| Sample | Average Density \pm Standard Deviation (g/ml) | | |
|------------------|---|-----------------|-------------------|
| | Pipette | Grad. Cyl. | Burette |
| <i>Coke</i> | 1.038 \pm 0.002 | 1.02 \pm 0.01 | 1.035 \pm 0.006 |
| <i>Diet Coke</i> | 0.997 \pm 0.007 | 0.98 \pm 0.02 | 0.995 \pm 0.004 |

(NOTE: Numbers listed in the Table above should be used as an example ONLY. They do not represent a real set of data. Note the use of the significant figures.)

Part 2. "Does the size of the sample affect the density?"

Sample: Diet Coke

Volume assigned: 20 ml

| | Trial 1 | Trial 2 |
|------------------------------------|---------|---------|
| Mass of Empty Beaker (g) | 29.0454 | 32.4122 |
| Mass of Beaker + Liquid Sample (g) | 48.9260 | 52.3385 |
| Mass of Liquid Sample (g) | 19.8806 | 19.9263 |
| Volume Delivered (ml) | 20 | 20 |

(NOTE: Numbers listed in the Table above should be used as an example ONLY. They do not represent a real set of data.)

| Sample | Density (g/ml) determined from the plot |
|------------------|---|
| <i>Coke</i> | 1.033 |
| <i>Diet Coke</i> | 0.995 |

Name Mark Prior

Date 09/14/05

TA Carlos Zambrano

Section # 52420

Discussion Questions:

Answer Discussion Questions listed at the end of the Procedure for Experiment 1.

(Example: The average densities and their corresponding standard deviations calculated from the class combined data for both Coke and Diet Coke are listed in the Table above.

Briefly, the following sets of data were obtained

for the Coke sample:

1.038 ± 0.002 (pipette)

1.02 ± 0.01 (grad. cylinder)

1.035 ± 0.006 (burette)

for the Diet Coke sample:

0.997 ± 0.007 (pipette)

0.98 ± 0.02 (grad. cylinder)

0.995 ± 0.004 (burette)

I would expect the density of Coke to be higher than the density of Diet Coke due to its higher sugar content...)

NOTE: *Attach a print-out of your data table and three graphs to your report.*