



**CH204 Experiment 1**  
**Dr. Brian Anderson**  
**Spring 2008**

**Are the Densities of  
Coke and Diet Coke  
Different?**

  $d = \frac{m}{V}$  

---

---

---

---

---

---

---

---

**Today**

- **Error in experimental data**
  - Random
  - Systematic
  - Gross
- **Accuracy and precision**
  - Accuracy - how close your final answer is to the correct one
  - Precision - how close your data points are to each other

---

---

---

---

---

---

---

---

**More Today**

- **Significant digits**
  - Count 'em!
  - Add and subtract 'em!
  - Multiply and divide 'em!
- **Standard deviation**
  - A statistical measure of random error
- **Quick look at Experiment 1**

---

---

---

---

---

---

---

---

**How much does a quarter weigh?**



---

---

---

---

---

---

---

---

**Weigh a few more...**

5.7128	5.6947
5.7085	5.6907
5.6106	5.6339
5.6009	5.7205
5.6466	5.7195

Now what does a quarter weigh?

**Average = 5.67387 grams**

---

---

---

---

---

---

---

---

**There's got to be a better way!**

**There is. Calculate the average  $\pm$  standard deviation**

**Standard deviation is a measure of the scatter in the data - and it also tells us how many decimal places we should keep in our average.**

---

---

---

---

---

---





---

---

## What is Standard Deviation?

It's a calculation based on a set of data points that tells us how widely the data points are scattered around the average.

$$\sigma = \left[ \frac{\sum (x_i - \bar{x})^2}{(n-1)} \right]^{1/2}$$



---

---

---

---

---

---

---

---

## Calculating Standard Deviation

This is 2008. Don't calculate it by hand. Use a built-in calculator function or use Excel.

Let's head to Excel right now and see how E-Z this is.

---

---

---

---

---

---

---

---

## So now what does a quarter weigh?

5.67287 ± 0.046377 g

Variability occurs HERE

Report the final answer only up to the first uncertain digit - which is the same decimal place where the standard deviation falls.

---

---

---

---

---

---

---

---

**This is really important, folks!**

$$5.67287 \pm 0.046377 \text{ g}$$

Round the standard deviation to ONE significant digit: 0.05

And report the average only up to that decimal place: 5.67

---

---

---

---

---

---

---

---

**Variability (random error) limits your answer**

$$5.67287 \pm 0.046377 \text{ g}$$

should be reported as

$$5.67 \pm 0.05 \text{ g}$$

And that's what a quarter weighs!

---

---

---

---

---

---

---

---

**So what does a quarter weigh?**

$$5.67 \pm 0.05 \text{ g}$$

68% of all quarters should weigh between 5.62 and 5.72 grams.

5.7128	5.6947
5.7085	5.6907
5.6106	5.6339
5.6009	5.7205
5.6466	5.7195

---

---

---

---

---

---

---

---

**Variability limits  
significant digits**

**There were five significant digits in the mass of each quarter, but only three significant digits in the final result.**

**The last two digits are insignificant because they are less than the variability in the measurement.**

**“Variability in the measurement” = random error.**

---

---

---

---

---

---

---

---

**Ways of Determining  
Random Error**

**For a single reading:**

- Precision of the equipment**
- Tolerance of the glassware**

**For many readings:**

**Statistics**

**That’s what we’re gonna do in lab today.**

---

---

---

---

---

---

---

---

**Std. Dev. = random error**

**We’ll use this same procedure to determine random error and significant digits in Experiment 1.**

**And speaking of Experiment 1...**

---

---

---

---

---

---

---

---

## Equipment


Graduated cylinder 0.1 ml

Volumetric pipette 0.01 ml

Burette 0.01 ml

Analytical balance

0.0001 grams!



---

---

---

---

---

---

---

---

## Two-Part Lab

### Part One:

- Measure the mass of 5 mL of sample using the analytical balance and three different types of glassware (pipette, burette, and graduated cylinder).
- Calculate density. Total of six data points.
- Enter your results into the spreadsheet on the computer nearest the printer, and use all the class data in your report.

---

---

---

---

---

---

---

---

## Two-Part Lab

### Part Two:

- Dispense your assigned volume using a burette, and measure the mass of the sample on the analytical balance.
- Do NOT calculate density.
- Enter your mass and volume measurements into the spreadsheet on the computer nearest the door, and use all the class data in your report.

---

---

---

---

---

---

---

---

## Important!

You will need all three graphs:

**Part One:**

**1 - Density chart and graph comparing different methods (includes average and standard deviation for each method).**

**Part Two:**

**2 - Mass vs volume graph for Coke**

**3 - Mass vs volume graph for Diet Coke**

---

---

---

---

---

---

---

---

## Handling bad data

**If you know it's bad - because you know something went wrong, or because the number is physically impossible - you can discard it.**

**If you don't like it because it's widely scattered, you can't just toss it, you have to apply the Q-test (see the appendix of the lab manual).**

---

---

---

---

---

---

---

---

## Interpolation

**In order to calculate the density of water at the same temperature as your Coke or Diet Coke sample, you will have to **interpolate** between the density values in the table on page 8 of the lab manual.**

---

---

---

---

---

---

---

---

**Final comments**

**When entering data, type with your fingers, not with your thumbs.**

**Beakers are not volumetric!**

**Remember to rinse your burette and put it away.**

---

---

---

---

---

---

---

---

**Final final comments**

**Next week: Final Exam, Part 1.**

**There is a sample quiz on the web site [Freebies](#) page.**

**Also play with the spreadsheet on the Freebies Page.**

***Bring a calculator!***

---

---

---

---

---

---

---

---